

Passwords required changing time

Password for **DASH 2000, 3000, 4000, 5000** and **Solar Monitors**:

[Enter day and month from monitor screen with leading zeros. \(Example: July 4 = 0407\).](#)

Reference Service Manuals: **DASH 3000, 4000 and 5000**; 2000966-542, **Solar 9500**; 2001085-078, **Dash 2000**; 2005873-010, **Solar 8000m and i**; 2026265-075, **Solar 8000**; 414993-001.

Passwords for DASH 2500:

[No Password required changing time.](#)

Passwords for B650 and B850 series monitors:

[User: biomed, Password: Change Me](#)

Reference Service Manuals: **B650**; M1165897-006, **B850**; 2040384-004.

Passwords for B20, B30 and B40 series monitors:

[No Password required changing time.](#)

Passwords for V100 and Pro Series monitors:

[No Password required changing time](#)

Passwords for CIC Version 4:

MS-DOS Password: [MMS_COM](#)

Reference Service Manuals: **CIC V4 (Bedrock)** Service Manual; 2001099-156.

Passwords for CIC Version 5:

[User: biomed, Password: Change Me](#)

Reference Service Manuals: **CIC (MP 100)** Service Manual; 2026419-130, **CIC (Bedrock)** Service Manual:
2026419-131
